




ÉCOLE POLYTECHNIQUE
FÉDÉRALE DE LAUSANNE



LABORATORY FOR BIO AND
NANO INSTRUMENTATION

USER GUIDE

LBNI AFM Controller Software

USB-7856R-OEM

Charlène Brillard

05.10.2018 – Beta Version 2.0.8

This Guide is a beta version

Not finished

Will be improved soon

For any question, contact me at familynamefirstname@gmail.com

ENVIRONMENT: LabVIEW 32bits



The Code has been developed with Labview 2016 (32 bits).

The Code can run on several targets

USB-7856R OEM

USB-7856R

PXIe-7975R

PXI-7851R

PXI-7841R

Full code


Feedback

Scan

But the easiest to interface is the
USB-7856R OEM

When debugging, check your signals
with an external oscilloscope

Open the Project LBNIController.lvproj



The software with his modules.

Menu.vi is the VI to run at first and will offer an interface to open the other modules.

Functional Global Variable.

If you want to change parameters at the initialization like the Calibration.

The target with the FPGA Code

A bitfile containing the FPGA Code optimization is loaded on the target by the host.



When a VI is locked: replace the name of the VI and run

Host: the top-level code on the computer

FPGA: The low-level code on Limited Physical Resource for Real-Time operations

MENU


From the project, open FrontPanel\Menu.vi and run the vi.



Navigate through the menu to select the modules you want to work with

Click on your choices, the modules windows will be opened

MENU – 10 MHz DAC/ADC



Change the clock signal value and switch from the Analog Input/Output to Custom ADC/DAC
To get a higher Sampling Rate for the Feedback and Excitation Signals

Available when there is a DAC/ADC attached to the respective DIO lines.

ZCONTROLLER

In the Submenu « Modules »

Feedback modes

- Contact
- Off-Resonance Tapping
- Auxiliary: user configurable

Enable the PID/ Disable the PID = enabling the manual motion

/!\ Disabling the PID doesn't retract the piezo. It freezes the position or enable to manually change the Z piezo position /!

Press on these buttons to move the motor up / down
The longer you press, the most steps


Time between two steps : 1ms is standard (min)
It makes sense to change it if you want to change the position of the piezo or approach manually

Start/Stop the Auto-Approach to find the surface

Number of steps done by the stepper motor during the auto-approach between two enabling/disabling of the Feedback. Multimode: 1 step = 200-300nm

To set the speed of the auto-approach, **change the I Gain** and/or the setpoint.

Blue led: Auto-Approach is running
Green led: Surface found



Number of steps moved by the Stepper motor

Fully retracts the piezo.

Press 4s and the stepper motor will move of the number of steps in «Approach»

PID on: Indicates the position of the Z piezo
PID off: Control the position of the Z piezo

Dragging the black cursor in the blue bar, clicking in the bar, or replacing the digital indicator by the requested value move manually the Z piezo to the requested position at a speed set in «Other Config»

Set the parameters of the PI Controller and the setpoint in V

Move the tip up relative to the current position by the Lift height

Tip: Before starting the auto-approach, check how fast the piezo extends with the current gain and setpoint by enabling the feedback, for a soft approach

ZCONTROLLER

- Other config


Select the type of system used to approach the tip and the sample. The left panel of Z Controller will adapt depending on the needed parameters

- Stepper Motor (Multimode : A, B, C, D digital signals)
- DC Motor
- Stepper Motor Driver (Enable, Direction, Steps)
- Stick-Slip stage (+ open the complementary module)

Set the direction for the up/down motor movement (the config by default is pretty standard)


Auto-Approach configuration:

- Do you want a beep sound when the surfaces is found
- Do you want to stay landed on the surface or withdrawn
- Do you want to hack the auto-approach (slow)



Relative enable: Engage on a percent decrease of the setpoint (useful for drifty signals, can find application for tapping mode)

SCAN



Changing the pixel/line restarts the scan from zero

If activated, save the picture at the end of each scan
Base Folder name (should be an existing path) and finish by \




Create a folder for the experiment at every name update.
Placeholders <date>_exp<#> can be removed or completed

The screenshot shows the 'Scan Control' window. It has a title bar with standard window controls. The main area is divided into several sections: 'Scan Parameters' with fields for Width (50um), Height (50um), X Offset (0m), Y Offset (0m), Angle (4.00°), Pixels (256), Lines (256), Scan Rate (1Hz), and Speed (100um/s); a 'Frame' section with buttons for Frame down, Frame up, Scan the same line, Alternate direction, Pause, Lock, Stop, and Change speed; a 'Saving' section with 'Auto save' checked, a 'Now' button, a folder path field (C:\), a file name field (<date>_exp<#>), and a number field (0); and a status bar at the bottom showing 'Scan.Disable finished'.



- Frame down (at every beginning of image)
- Frame up (at every beginning of image)
- Scan the same line in case of false and make a picture of it
- Alternate the direction of the scan
- Pause the scan at the current tip position.
To continue the scan, press Pause again.
- Lock the value to the same of the top one
- Stop the scan and move the tip to zero position
- Change the speed display parameter from Hz to s(/line)
- Save the current picture (referring to what is displayed in the monitor)
- Browse the Folder
- # of the experiment (replace the placeholder)
- Tells you what's happening, it is busy when not written « ...finished »

MONITOR


Open as many as needed from «Graphs» by clicking several times on «Monitor»



MONITOR



MONITOR



Update the graph with the current parameters

«Rounded» triangle is the default choice for the fast axis. An abrupt change of direction in the scan (triangular waveform) creates ripples particularly visible at high scan speed. Rounding helps to reach higher speeds.

To ignore for the current development


Displays the waveform with the hysteresis correction on top

Apply the lateral resonance compensation from the filter

Indicate if the loaded waveform contains the correction

Send the waveform displayed in the graph to the FPGA

Generate 3 types of waveforms




Tip: The higher the scan speed, the more rounding is needed to avoid the ripples. Rounding also reduces the scan range as the voltage ramp applied to the piezo is modified.


SIGNALS

Ctrl+click to select the channels to save
(Forward and/or backward)

Ctrl+click to select channels to display
Resets the channels to save



Change the text by clicking once in a line



Choice of window panel: Selection/Overview/Assignment

Change the unit and calibration of the displayed and saved data

In **Overview**, monitor some signals (under development)

Signals					
Selection		Overview		Assignment	
Ch	Name	Value	Unit	Cal	Offset
AI0	Current	630.5m	m	1	0
AI1	AI1	-1.734	V	1	0
AI2	AI2	2.484	V	1	0
AI3	AI3	-10	V	1	0
AI4	AI4	-10	V	1	0
AI5	AI5	9.239	V	1	0
AI6	AI6	7.294	V	1	0
AI7	AI7	-10	V	1	0
AO0	X	4.084u	m	1	0
AO1	Y	-6.337u	m	1	0
AO2	Z	2.265u	m	1	0
AO3	AO3	0	V	1	0
AO4	Bias	0	V	1	0
AO5	DC Motor	0	V	1	0
AO6	AO6	0	V	1	0
AO7	AO7	0	V	1	0

RAMP

Ctrl+click to select channels to display and save

There will always be a Backward

Can enable and disable the Forward movement for Pull

Choose the channel to ramp: Z (Provision for Bias ramp for example)

Reenable the feedback after the ramp

Calculate the deflection sensitivity (in nm/V) by calculating the slope between A and B cursors on the error channel

Has to be an existing path (Read from the path updated from Scan.vi)

If empty, replaced by Ramp_date_x.txt

To play with bigger number of points, reduce the recordClockDec (careful with labview memory, not optimized yet)

Save the current displayed force curved

Start Ramping with the current parameters

Not useful for now

Number of ramps which will be done

Stops the ramp

Ramping state

Absolute Height: read positions from piezo

Process Variable: signal used for the feedback (deflection setpoint for example)

External: not tested yet

Relative height: start = the distance withdrawn from current position. distance = the ramp distance (from start)

- = from the surface
+ = towards the surface

Number of points
If red, too high for the record clock (advices given in context help)

Speed of Fwd movement

Lock the Bwd to Fwd speed
Delay after feedback is reenabled (useful for several ramps)

Blue = Forward \\ Red = Backward (Sensitivity calibration: 1/slope)

EXCITATION

Changing the Frequency, the amplitude, the DC Offset (disables the excitation, withdraws the cantilever) loads the new parameter (and reenables the excitation and the Feedback). The Enable thick-mark is changing of status and text during this process.

ORT Piezo : Signal added to the Z available in channel AO2. For example: $AO2 = Z + \text{Amplitude} \times \sin(Ft) + \text{DC Offset}$
 PORT Photothermal: Signal sent to the Drive in channel AO3. $AO3 = \text{Amplitude} \times \sin(Ft) + \text{DC Offset}$

Choice between displaying the Excitation signal or the Hydrodynamic Background after Correction

Sync point is where on the interaction curve the feedback is done.




Before engaging, the *SyncPoint* has to be on and set to a sensible value, from the background shape (automatically set at the minimum in case of PORT, can be refined. In case of ORT, guessed from experience.)

Average width around the *SyncPoint* in points where the feedback is done.

Correct the background at every *end of frame*, *end of line* or *never*.

Number of background curves taken to calculate the background and subtract it to the feedback signal. The higher, the longer it takes to calculate

Distance withdrawn from the surface on correcting background



EXCITATION

Photothermal


Enable the excitation in the drive channel (AO3)

Frequency of the signal


The value entered is coerced to the value really applied

Limited by the sampling rate of the analogs channels of the FPGA: With 1MS/s, $F=30\text{kHz} \Rightarrow 33$ datapoints

Frequency of the signal




Excitation Signal





Cuts of the excitation by a certain pourcentage to reduce the duration of the pulse, particularly useful when exciting the cantilever thermally or photothermally)

Shits the excitation pulse in the memory for visualisation of the interaction curve



EXCITATION





Height




Normal limits of the piezo :
 Lower = -Piezo range / 2 (-2.5um here)
 Upper = +Piezo range / 2 (+2.5um here)

To avoid clipping and damaging the piezo,
 Limits are defined in ORT:
 Lower = -Piezo range / 2 + ORT Amplitude
 Upper = +Piezo range / 2 - ORT Amplitude

ORT FORCE CURVES



Filtering



OSCILLOSCOPE

